

SciCan OPTIM

One-Step Cleaner-Disinfectants



DEMONSTRATED VIRUCIDAL EFFICACY EXTENDS TO NEW & EMERGING PATHOGENS INCLUDING COVID-19

The global COVID-19 pandemic has affected supply chains in many different industries. In response to the record-breaking demand for OPTIM disinfectants, we have ramped up production. We continue to produce and ship OPTIM Wipes and Solution to dental dealers every week.

The global supply of raw materials used to manufacture pre-saturated wipe products has been substantially constrained as a result of the pandemic. This means you might be experiencing longer than usual lead times on these items. If you cannot obtain OPTIM Wipes, consider using OPTIM Liquid. After all, our Wipes and Liquid use the same highly effective formula.

OPTIM has demonstrated effectiveness against viruses similar to the COVID-19 coronavirus on hard non-porous surfaces and non-invasive medical devices. This product can therefore be used against COVID-19 when used in accordance with the directions for use against Poliovirus on hard, non-porous surfaces, and non-invasive medical devices.



Using OPTIM Liquid

1. Spray SciCan dry wipe with OPTIM Solution.
(Dry Wipes - Item Number: SCI33R18 / SCI33R18XL)
2. If surface is not heavily soiled, clean and disinfect surface with one wipe and allow surface to remain wet for the contact time indicated on the label.
3. If surface is visibly soiled, clean surface with one wipe and discard. Use a second wipe to disinfect and allow surface to remain wet for the contact time indicated on the label.




OPTIM F.A.Q.



Is it OK to place disinfectant solution in a container with 4X4 gauze for use on dental equipment?

Use of pre-soaked gauze or paper towel is not recommended. Gauze is naturally brown in colour. The only way to make it white is to heavily bleach it. Presoaking gauze in a disinfectant solution implies putting two competing chemistries head to head in a container which is never a good idea. All gauze tends to absorb too much product and can cause pooling on surfaces, which is inadvisable. Use only SciCan Dry Wipes Item Number SCI33R18 / SCI33R18XL.



What is OPTIM's contact time for Novel Coronavirus (COVID-19) and what does that mean?

The contact time is the length of time a disinfectant must be wet on a surface to achieve the desired efficacy result. OPTIM 1 has a 30-second virucidal contact time. That means that the surface must remain wet for only 30 seconds in order to inactivate viruses including the Novel Coronavirus (COVID-19).

Can OPTIM be used to clean surfaces?

Yes, OPTIM has been tested as a one-step cleaner disinfectant, meaning it cleans your surfaces as it disinfects. Please refer to the product label to learn how to use our product as a one-step cleaner-disinfectant.



Which surfaces can be cleaned and disinfected by OPTIM?

OPTIM was designed for use on hard, non-porous surfaces and non-invasive medical devices. This includes dental chairs, stools, lights and counter tops if used according to our instructions for use on the label. Please refer to the Instructions for Use on the product label.

OPTIM Ordering Information

Pre-Saturated Wipes

OPTEU1-W12 OPTIM 1, Ready-to-Use Wipes – 6" x 7" (12 x 160 sheets)

Liquid

OPTEU1-4X4 OPTIM 1, Ready-to-Use Liquid (4 x 4L)

Accessories

FSA01 Sprayer for 1 L Bottle

SCI33R18 SciCan Dry Wipes – 6" x 7" (18 x 160 sheets)

SCI33R18XL Large SciCan Dry Wipes – 10" x 10" (18 x 60 sheets)

Dent4You AG
Bahnhofstrasse 2
CH-9435 Heerbrugg



EC REP

Coltene/Whaledent GmbH+Co. KG
Raiffeisenstraße 30
DE-89129 Langenau