

SciCan STATIM G4

Cassette Autoclaves



**Sterilization and Dryness at speeds faster than
conventional chambered autoclaves**



Fast Cycle Times

STATIM 2000 G4 can sterilize unwrapped loads in under 7 minutes, and wrapped loads in under 11 minutes.

STATIM 5000 G4 can sterilize unwrapped loads in under 9 minutes, and wrapped loads in under 16 minutes.



Reliable

For over 25 years, STATIM's unmatched speed, reliability, and dependability have made it one of the world's best selling cassette autoclaves.



Capacity

STATIM G4 autoclaves not only offer fast sterilization times but also enough capacity to meet the needs of a busy practice.

STATIM 2000 G4 capacity allows for up to 4 pouches and the STATIM 5000 G4 allows for up to 10 pouches.



G4 Technology

SciCan's G4 equipment can connect and send information to your smart devices so you can have cycle data, instructions, and maintenance notifications at your fingertips.



Gentle

The fresh water feature guarantees fresh water with every cycle. State of the art steam technology is gentle on both solid and hollow instruments.



STATIM 2000 G4 Specifications

Cycle Name	Load	Sterilization Temperature	Sterilization Time	Total Time (mm:ss)*	
				Warm Start	Cold Start
SOLID UNWRAPPED (N)	1.0kg	134°C	3.5 min	6:45	9:15
HOLLOW UNWRAPPED (S)	1.0kg	134°C	3.5 min	8:05	11:45
HOLLOW UNWRAPPED (S)	1.0kg	134°C	18 min	22:35	26:15
HOLLOW WRAPPED (S)**	1.0kg	134°C	3.5 min	10:40	15:35
HOLLOW WRAPPED (S)	1.0kg	134°C	18 min	25:10	30:05
RUBBER/PLASTIC (S)	0.4kg	121°C	15 min	18:40	20:15
RUBBER/PLASTIC (S)	0.4kg	121°C	30 min	33:40	35:15
AIR DRYING ONLY			Not a sterilization cycle		

Unit size (L x W x H):
19.5" x 16.25" x 6" / 49.5 cm x 41.5 cm x 15 cm

Cassette internal dimensions:
11" x 7" x 1.4" / 28 cm x 18 cm x 3.5 cm

Minimum shelf height for proper leveling:
8.5" / 22 cm

Weight without water:
46 lbs / 21 kg

Power consumption:
220-240 V, 50/60Hz, 6A

* Drying not included
**Cycle for routine check (STATIM PCD test)

STATIM 5000 G4 Specifications

Cycle Name	Load	Sterilization Temperature	Sterilization Time	Total Time (mm:ss)*	
				Warm Start	Cold Start
SOLID UNWRAPPED (N)	1.5kg	134°C	3.5 min	8:45	13:15
HOLLOW UNWRAPPED (S)	1.5kg	134°C	3.5 min	10:50	17:30
HOLLOW UNWRAPPED (S)	1.5kg	134°C	18 min	25:20	32:00
HOLLOW WRAPPED (S)**	1.5kg	134°C	3.5 min	15:30	24:00
HOLLOW WRAPPED (S)	1.5kg	134°C	18 min	30:00	38:30
RUBBER/PLASTIC (S)	0.4kg	121°C	15 min	20:20	22:50
RUBBER/PLASTIC (S)	0.4kg	121°C	30 min	35:20	37:50
AIR DRYING ONLY			Not a sterilization cycle		

Unit size (L x W x H):
23.6" x 16.1" x 7.5" / 60 cm x 41 cm x 19 cm

Cassette internal dimensions:
15" x 7" x 3" / 38 cm x 18 cm x 7.5 cm

Minimum shelf height for proper leveling:
10" / 26 cm

Weight without water:
73 lbs / 33 kg

Power consumption:
220-240 V, 50/60Hz, 6A

* Drying not included
**Cycle for routine check (STATIM PCD test)



VistaPure

Water Filtration

Provides two grades of high-quality water for autoclaves, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers.

V3000



STATMAT

Protect your STATIM

Securely seat cassettes atop your STATIM with a safe and hygienic silicone mat.

01-112867S - For STATIM 2000

01-112868S - For STATIM 5000



STAT-DRI Plus

Drying Agent

Enhance drying in your STATIM Cassettes with the STAT-DRI Plus drying agents.

2OZPLUS - Stat-Dri Plus 2 oz with Sprayer

8OZPLUS - Stat-Dri Plus 8 oz with Sprayer

8OZPLUS - Stat-Dri Plus 8 oz with Cap (refill)

32OZPLUS - Stat-Dri Plus 32 oz with Cap (refill)

T473 - Spray Trigger Only for 32 oz Bottle

Dent4You AG
Bahnhofstrasse 2
CH-9435 Heerbrugg

Manufactured by:
SciCan Ltd.
1440 Don Mills Rd.
Toronto, ON, M3B 3P9 / Canada,
T +1 416-445-1600
TF +1 800-667-7733
customerservice@scican.com

VistaPure™ is manufactured by:
Crosstex International, Inc.
1733 W. Parkside Lane
Phoenix, AZ 85027 / USA
T +1 800-531-2221
VistaResearchGroup.com

www.scican.com

COLTENE