STATIM® StatMat



Keep your STATIM looking new with the StatMat Silicone Mat



Provide a safe and hygienic working area by securely seating the cassette on top of the STATIM.

Protects paintwork from accidental scuffs and scrapes, keeping your STAT/M looking new for longer.

Helps dissipate heat and protects surfaces from damage caused by heat transfer from hot metal surfaces.

StatMats are available for both STAT/M 2000 G4 (01-112867S) and STAT/M 5000 G4 (01-112868S) autoclaves.





Code	Description
01-112867S	Silicone Mat STAT/M 2000S/G4
01-112868S	Silicone Mat STATIM 5000S/G4
01-1107548	Silicone Mat STAT <i>IM</i> 7000

